



**NORTH AMERICAN IMAGING**

## ***Single Piece Integrated 9-inch X-ray imaging system***

This new imaging system, exclusively from NAI, is an integrated system equipped with a 9-inch image intensifier and a CCD camera.

### **System Features:**

- Single piece streamlined housing design
- High-resolution SD series image intensifier
- Noise reduction and motion detector image processing technology
- Full DSP digital image processing function (digital signal processing)
- Compact lens with optional automatically-controlled iris mechanism
- Distortion Correction Lens
- Eight-frame memory
- OSD (on-screen display) function



**Contact X-Ray Associates, LLC**  
**1-909-599-9150**



**X-Ray Associates, LLC**

482 West Arrow Highway, Unit L  
San Dimas, CA 91773-2941  
Phone 909-599-9150 Fax 909-599-9160  
<http://www.xrayllc.com>